



Hearing Systems – Part 3 of 3

<https://hearingsystems.github.io/>


TU Ilmenau – Audio Signal Processing & Audio Systems

January 24, 2024


Dr.-Ing. Tamas Harczos

•• Recap: Equal-loudness contours

Equal-loudness contours




Exemplary
Equal-loudness contours of a
hearing impaired (HI) individual



Based on: Stroebe & Alwan, 1995


•• Hearing aids: Wide dynamic range compression



Linear amplification alone:

- Increases amplitude always with the same gain.
- Not conform with loudness perception of hearing impaired. (Will often cause uncomfortably loud percept.)
- Can lead to more hearing impairment.


•• Hearing aids: Wide dynamic range compression



Compression:

- Example:
40 dB increase in the input leads to 20 dB increase in the output, i.e.,
→ **compression ratio: 40 dB : 20 dB = 2:1**
- Compression ratio can vary, but usually does not exceed 4:1 in the hearing aid field.


•• Hearing aids: Wide dynamic range compression



Compression:

- Low level gain: gain at ca. 50 dB SPL input
- High level gain: gain at ca. 90 dB SPL input


•• Hearing aids: Wide dynamic range compression



Second linear range:


- Can be useful to allow the listener to realize that the input is very loud.
- Should be verified to not extend into the threshold of pain. (See “limiter” on next slide.)
- Remember the normal hearing ear handles loud noises very similarly! Basilar membrane nonlinearity:

Reminder



Adapted from: [Johnstone et al., 1986](#)


•• Hearing aids: Wide dynamic range compression



Limiter:

- Ensures not to cross the threshold of pain.
- Ensures technical limits not to be exceeded.
(MPO: maximum power output of a hearing device. Limited for example by the power source, amplifier, and receiver.)
- Is equivalent to a compressor with $\infty:1$ compression rate.


•• Hearing aids: Wide dynamic range compression



Expansion:


- Usually employed in the low SPL range.
→ Less gain added to soft noises.
- Reduces inherent noise of the system (e.g. that of microphones or preamps).
- Reduces irrelevant background noise.
→ Increases overall SNR.
- Typical expansion ratios: 1:1.5 - 1:3.
(1:1.66 in the graph on the left.)

•• Hearing aids: Wide dynamic range compression



Gain plot:

- A shorter / compacter way to describe compression-expansion curves.




•• Hearing aids: Multi-channel WDRC

Multi-channel compression: Why?

- The dynamic range of hearing (normal or impaired), the amount of hearing loss, and also the dynamic range of speech [KR16] vary across frequency.

Multi-channel compression: How?

- Using multiple dynamic compressors, each fitted for a distinct frequency range.
- A computationally cheap, but high quality filter-bank with almost perfect reconstruction is needed. E.g. an FFT-based WOLA filter bank [RC80].



[KR16] M. Kirchberger & F. Russo, "[Dynamic Range Across Music Genres and the Perception of Dynamic Compression in Hearing-Impaired Listeners](#)," *Trends in Hearing*, vol. 20, no. 4, doi: 10.1177/2331216516630549, 2016.

[RC80] R. E. Crochiere, "[A weighted overlap-add method of short-time Fourier analysis/synthesis](#)," *IEEE Trans. Acoustics, Speech and Signal Proc.*, vol. 28, no. 1, doi: 10.1109/TASSP.1980.1163353, 1980.

•• Hearing aids: Multi-channel WDRC

Multi-channel compression


- The dynamic range of the hearing aid is compressed to match the dynamic range of speech

Multi-channel compression


- Using multiple dynamic range compression across distinct frequency ranges
- A computationally complex algorithm with almost perfect frequency selectivity (e.g., FFT-based WOLA filterbank)

[KR16] M. Kirchberger & F. K. ...
in Hearing-Impaired

[RC80] R. E. Crochiere, "Acoustics, Speech and Signal Processing, vol. 2, pp. 1301-1308, 1980."



s, and also the



Dynamic Compression
6630549, 2016.

"IEEE Trans.

•• Recap: Causes of hearing impairment



- Malformation
- Foreign body




- Eardrum scars / malformation / perforation
- Barotrauma (can be temporary)
- Otitis media or Eustachian tube dysfunction (fluid in the middle ear)
- Rigid joints between ossicles, otosclerosis



- Aging (presbycusis)
- Exposure to loud noise
- Trauma / Barotrauma / Head injury
- Genetic or hereditary hearing loss
- Viral infection of the hearing nerve
- Malformation of the inner ear or tumor
- Ototoxic drugs (antibiotics, some anti-inflammatory drugs, chemotherapy)



- Impacted ear wax
- Foreign body or fluid
- Inflammation or infection
- Malformation or tumor



Source: www.lobe.ca, 2021.

•• Cochlear implants: Early history

Alessandro Volta
50V DC ear-to-ear




Source: hearinghealthmatters.org, 2021

Lundberg (1950),
Djourno & Eyries (1957)

Direct electrical stimulation of the auditory nerve with distinguishable percepts in human.

House's first clinical
multi-channel CI



Source: [Clark](https://www.clark.com), 2013

Clarion implant
Advanced Bionics
(USA)



Source: [cochlear.com](https://www.cochlear.com), 2021

1790s

1868

1950-57

1961

1982


1985

1996

2001


2006

Source: abebooks.com, 2021



Rudolf Brenner's
AC stimulation

Source: [Elsenberg](https://www.elsenberg.com), 2005



William F. House's
single-channel CI




Nucleus implant
Cochlear Ltd. (AU)

Combi40+ implant
MED-EL GmbH (AT)


Venus CS-10A implant
Nurotron BioTech (CN)

Commercialization →

•• Cochlear implants



- Cochlear implants




Source: [cochlear.com](https://www.cochlear.com), 2021

- Cochlear implants



Source: [cochlear.com](https://www.cochlear.com), 2021 (unprocessed sound)

•• Cochlear implants: The electrode array




Source: earsurgery.org, 2021

Misplacement of the cochlear implant (CI) electrode array can happen. See e.g.:


G. Widmann, D. Dejaco, A. Luger et al., "[Pre- and post-operative imaging of cochlear implants: a pictorial review](#)," *Insights Imaging*, vol. 11, no. 93, doi: 10.1186/s13244-020-00902-6, 2020.

•• Cochlear implants: Electrical stimulation

- Electrical stimulation must be **charge-balanced**!
- Otherwise :
 - electrolysis happens (fluid*-filled cochlea);
 - chemical balance gets corrupted;
 - cells die in the cochlea.
- Most widespread is the **biphasic** (negative first) **charge-balanced pulsed** type stimulation.
- Other types of multiphasic pulse shapes:
 - triphasic,
 - precision triphasic.



Source: [Tognola et al., 2005](#)




Source: [Abdelhamed, 2019](#)

*: Endolymph is a physiological fluid that fills the inner ear's labyrinth


•• Cochlear implants: Current spread (or *spread of excitation*)

- Cochlea is filled with endolymph.
- An electric pulse (independent of its physical properties) from a single electrode cannot be focused.
- Current will spread out along the cochlea.
- Current will also spread (to some extent) through the bony wall of the cochlea towards neighbouring turns.
- This means, we will always stimulate a group of auditory nerve fibres and not just the ones in the direct proximity of each CI electrode contact.




Source: [Jeschke & Moser, 2015](#)

•• Cochlear implants: CIS signal processing strategy




Continuous Interleaved Sampling (CIS) signal processing strategy



B. S. Wilson, C. C. Finley, D. T. Lawson, R. D. Wolford, and M. Zerbi, "[Design and evaluation of a continuous interleaved sampling \(CIS\) processing strategy for multichannel cochlear implants](#)," *J. Rehabil. Res. Dev.*, vol. 30, no. 1, pp. 110–116, 1993.

•• Cochlear implants: ACE signal processing strategy




Advanced Combination Encoder
(ACE) signal processing strategy


The most widespread CI strategy today.

A. E. Vandali, L. A. Whitford, K. L. Plant, and G. M. Clark,
 "Speech Perception as a Function of Electrical
 Stimulation Rate: Using the Nucleus 24 Cochlear Implant
 System," *Ear and Hearing*, vol. 21, no. 6, pp. 608-624,
 doi: 10.1097/00003446-200012000-00008, 2000.


•• Cochlear implants: Stimulation patterns (cochleagram vs. electrogram)



•• Cochlear implants: Other signal processing strategies



ACE: advanced combination encoder, AIS: asynchronous interleaved stimulation, APS: advanced peak selection, CA: compressed analog, CI: cochlear implant, CIS: continuous interleaved sampling, F0-ACE: ACE incorporating F0 cues, FSP: fine structure processing, FS4p: FSP with parallel stimulation and four channels featuring channel-specific sampling sequences, HDCIS: high-definition CIS, HiRes: high resolution (CIS), MP3000: ACE incorporating simultaneous masking effects, MPIS: main peak interleaved sampling, NVJCIS: novoicing jitter CIS, PES: pitch excited sampler, PLCIS: phase-locked CIS, PPS: paired pulsatile sampler, QPS: quadruple pulsatile sampler, SAS: simultaneous analog stimulation, SineEx: sinusoid extraction strategy, SpecRes: spectral resolution strategy, SPEAK: spectral peak strategy



•• Hearing aids (HA) vs. Cochlear implants (CI)

- **System size**

BTE processor is about the same size for both, but CI has more system parts (RF link, implanted stimulator, implanted electrode array).

- **Power consumption**

Signal processing requires about the same power for both, but CIs need extra current for electrical stimuli and RF link.

- **Core signal processing**

Both do spectral decomposition via an analysis filter-bank, but while HAs need to synthesize audio on the output, CIs calculate parameter for electrical stimuli.

- **Additional algorithms / features**

Very similar: directionality, dynamic range compression, noise reduction, binaural algorithms, but CIs do not need feedback canceller. Both can be streamed to wirelessly and be remote controlled.

- **Current issues**

Both: speech in noise enhancement, power consumption and size, connectivity, user-friendly handling.
CIs: *reverb*, *pitch discrimination* (12/16/22/24 channels, does not matter ← *current spread*).


•• Cochlear implants: Ways to improve pitch sensitivity

- Train cochlear implant users

Use the brain's plasticity. Needed anyway in CI therapy, but has its limits.

- Enhance pitch changes

Smart preprocessing can increase speech melody changes in real-time ([PREX project](#), [Fraunhofer IDMT](#)).



F. Kuhnke, L. Jung, and T. Harczos, "[Compensating for impaired prosody perception in cochlear implant recipients: A novel approach using speech preprocessing](#)," in *Proc. 5th Int. Symp. Audit. Audiol. Res. (ISAAR 2015): Individual hearing loss - Characterization, modelling, compensation strategies*, The Danavox Jubilee Foundation, pp. 309–316, Nyborg, Denmark, 2016.

F. Kuhnke, "Implementation and evaluation of a real-time pitch range extension algorithm for cochlear implants," M.Sc. Thesis (Technical University Ilmenau, Department of Electrical Engineering and Information Technology, Institute of Media Technology), 66 pages, 2014.

•• Cochlear implants: Ways to improve pitch sensitivity

- Train cochlear implant users

Use the brain's plasticity. Needed anyway in CI therapy, but has its limits.

- Enhance pitch changes


Smart preprocessing can increase speech melody changes in real-time ([PREX project](#), [Fraunhofer IDMT](#)).




- Improve temporal pitch cues

Better stimulation strategies ([SAM project](#), [Dissertation](#), [Fraunhofer IDMT](#))

•• Pitch perception




•• Cochlear implants: SAM signal processing strategy



Stimulation based on
Auditory Modeling
(SAM) signal
processing strategy

T. Harczos, "[Cochlear implant electrode stimulation strategy based on a human auditory model](#)," PhD Dissertation, urn: urn:nbn:de:gbv:ilm1-2015000204, 2015.

•• Cochlear implants: Stimulation patterns (ACE vs. SAM)



•• Cochlear implants: Ways to improve pitch sensitivity

- Train cochlear implant users

Use the brain's plasticity. Needed anyway in CI therapy, but has its limits.

- Enhance pitch changes

Smart preprocessing can increase speech melody changes in real-time ([PREX project](#), [Fraunhofer IDMT](#)).



- Improve temporal pitch cues

Better stimulation strategies ([SAM project](#), [Dissertation](#), [Fraunhofer IDMT](#))



•• Cochlear implants: Ways to improve pitch sensitivity

- Train cochlear implant users

Use the brain's plasticity. Needed anyway in CI therapy, but has its limits.

- Enhance pitch changes

Smart preprocessing can increase speech melody changes in real-time ([PREX project](#), [Fraunhofer IDMT](#)).



- Improve temporal pitch cues

Better stimulation strategies ([SAM project](#), [Dissertation](#), [Fraunhofer IDMT](#))




- Improve place pitch cues (reducing current spread)


- Better electrode design (“modiolus-hugging” or “adhering” electrodes),
- Current steering (parallel stimulation over multiple electrode contact points)
- Optogenetic stimulation (optical cochlear implants).

•• Optogenetic cochlear implants

Electrical stimulation



Optical stimulation



Source: Jenschke & Moser, 2015


Challenges of optogenetic cochlea stimulation:

- **Light-sensitive neurons:** efficient, reliable and safe opsin* expression via *genetic manipulation* of the spiral ganglion neurons.
- **Stimulation technology** must be efficient, reliable and safe.

*opsin = light-gated ion channel

•• Optogenetics


(C) Ninghui Shi - Wikipedia



Take the gene for the light-sensitive protein (ChR) from *Chlamydomonas reinhardtii*


Promoter

Gene encoding opsin



Build **genetic construct**
(including promoter and gene encoding opsin)

(C) Jazziw - Wikipedia




Insert genetic construct into **virus** (e. g. AAV)



Transfect animals
(injection into nervous system)
or use transgenic animals




Insert **optical fiber** or **μLED-array** and apply light of opsin-specific wavelength




The **neurons** expressing the opsin can be **activated** with high temporal and spatial resolution

•• Optogenetics




Source: be.mit.edu, 2021

•• Optogenetic cochlear implants: Stimulation technology (optrodes)




•• Optogenetic cochlear implants: The Team




•• Optogenetic cochlear implants


Electrical stimulation



Optical stimulation



Research hardware




Challenges of optogenetic cochlea stimulation:

- **Light-sensitive neurons:** efficient, reliable and safe opsin* expression via *genetic manipulation* of the spiral ganglion neurons.
- **Stimulation technology** must be efficient, reliable and safe.

*opsin = light-gated ion channel

- Optogenetic cochlear implants: Demo with the research hardware



Music: Josh Armistead - Lazer Beams

Licence: Attribution 4.0 International (CC BY 4.0)


•• Hearing aid market

2022:

1 million CI systems sold worldwide until 2022 [FZ22]

vs.

20.25 million hearing aids sold in Europe during 2022 [EH22]



[FZ22] F.-G. Zeng, "[Celebrating the one millionth cochlear implant](#)," *JASA Express Lett.*, vol. 2, no. 077201, doi: 10.1121/10.0012825, 2022.

[EH22] European Hearing Instrument Manufacturers Association, "[Hearing aid sales](#)," url: <https://www.ehima.com/about-ehima/hearing-aid-sales/>, 2022.

•• Future of hearing systems

Hearing aids:

- convergence towards earbuds (similar form factor and features like noise-cancellation);
- more biometric sensors (like PPG) for health-monitoring;
- neuro-steering (EEG based attention decoding and filtering);
- better wireless connectivity (MFi, ASHA, LE Audio, Auracast) while keeping power consumption ultra low;
- remote fitting → self-fitting → remotely supervised AI-based auto-fitting;
- OTC (over-the-counter) hearing aids.

Cochlear implants:

- fully implantable devices (maintenance-free for 6-15 years);
- optical cochlear implants (with the help of optogenetics).

•• Future of hearing systems

Hearing aids:

- conventional
- more powerful
- neural
- better
- while
- remote
- OTC

Cochlear

- fully
- optical



AI's take on future of hearing aids

- Generated by
Webjourney.ai


- Prompt:
"close-up of a
future hearing aid
on a human ear,"



•• PrimeHA: a platform to develop and test hearing aid algorithms


Main features:

- ARM Cortex M4f @ 64 MHz, 256 kB RAM, 1 MB NVM.
- Current consumption $\sim 3\text{-}10$ mA with BLE 5.1 radio.
- Audio input: 2x PDM MEMS microphones.
- Extras: motion sensor and μSD -card slot on board.
- Made of commercial off-the-shelf components.
- Fully C-programmable under free-for-all license.



Current use:

- Firmware framework development.
- Testing and optimizing audio signal processing algorithms like feedback canceller, multi-band compressor, noise reduction.
- Learning lessons on acoustically engineered mechanical design.



[AH23] A. Hintermaier, I. Pieper, and T. Harczos,
"PrimeHA: not a hearing aid," VCCA, 2023.



Thank you very much!
Questions?

Dr.-Ing. Tamas Harczos
tamas.harczos@audifon.com

audifon GmbH & Co. KG