



# Hearing Systems – Part 2 of 3

<https://hearingsystems.github.io/>

TU Ilmenau – Audio Signal Processing & Audio Systems

---

January 23, 2024

Dr.-Ing. Tamas Harczos

## •• Introduction

- audifon:
  - audifon GmbH & Co. KG develops, produces and sells hearing systems;
  - headquarter located in Kölleda, Thüringen (also place of production);
  - R&D offices in Köln (Electronics) and Ilmenau (Embedded Software);
  - about 220 employees.
- KIND:
  - KIND GmbH & Co. KG is specialized in hearing acoustics and optics;
  - headquarters located in Großburgwedel near Hanover;
  - around 750 specialist stores and over 3,000 employees;
  - KIND and audifon are both owned by the family Kind.




## •• Introduction

- Reminder:
  - Three lectures and one seminar.
  - Useful links:
    - <https://hearingsystems.github.io/> (presentation slides and recordings),
    - <https://moodle.tu-ilmenau.de/course/view.php?id=125>.
- Topics we will cover today:
  - *Recap of the basics*: **units**, human hearing, **tonotopy**.
  - A few more words on **inner and outer hair cells** and **equal loudness contours**.
  - Properties and characterization of hearing impairment: **pure tone audiogram** and **speech in noise tests**.
  - **Hearing aids**: early history, modern day form factors, features and requirements.

## •• Recap of the basics: scales and units

- Decibel (**dB**) alone is not a unit, it's a scale. It describes a relationship between two levels.
- The unit **dB FS** (Decibel full scale) relates to the maximum output of the given system.
  - 0 dB FS means maximum signal amplitude for a given system.
  - 1/100, 1/10, 1/2 of maximum amplitude: -40 dB FS, -20 dB FS, ca. -6 dB FS, respectively.
  - Conversion:  $dB = 20 \log_{10}(mag) \leftrightarrow mag = 10^{(dB/20)}$ .
  - A dB FS value does not necessarily relate to the loudness of a signal.
- The root mean square (RMS) of a signal is related to the energy residing in it.
  - The RMS of a recorded audio signal correlates with its sound pressure level.
- The unit **dB SPL** relates to the reference sound pressure  $p_0 = 20 \mu\text{Pa}$  in air.
  - 20  $\mu\text{Pa}$  is often considered as the threshold of human hearing.
  - Distance (source to receiver) is important, 1 m is frequently used as standard distance for dB SPL measurements.
  - Doubling the distance reduces the SPL by about 6 dB (in free field).

- Recap of the basics: human hearing




Source: [youtube.com/@Signia-hearing](https://www.youtube.com/@Signia-hearing), 2024.

## •• Recap of the basics: tonotopy





Source: [HHMI, 2000](#).

**Tonotopy** is the spatial arrangement of where sounds of different frequencies are processed in the auditory nervous system.



## •• A few more words on tonotopy



### Frequency scale of the ear :


- The human basilar membrane (BM) is approx. 32-35 mm long. Along the BM ca. 3500 IHCs and 12000 OHCs are distributed.
- Regardless of whether we take Bark or Mel scale, it is important to remember that place pitch is not distributed linearly along the cochlea;
  - 210±10 Hz ( $\Delta=0.2$  Bark): good discrimination, but
  - 4410±10 Hz ( $\Delta=0.03$  Bark): bad discrimination.




- Pitch difference of 1 Bark roughly equals to 1.3 mm on the BM; you only have around 150 IHCs along that cochlear section.

## •• How hearing works: Inner hair cells

### Organ of Corti



Adapted from: [Baumgarte, 2000](#)



Source: [Brown & Santos-Sacchi, 2014](#)


### Inner hair cells:

- do mechanoelectrical transduction (converting motion to electricity);
- are the sensory cells in our auditory system.



## •• How hearing works: Outer hair cells

### Organ of Corti




Adapted from: [Baumgarte, 2000](#)

### OHC motility



Source: [J. Ashmore for BBC 'Ear we go', 1987](#)

### BM nonlinearity




### Outer hair cells (OHCs):

- actively amplify/dampen the movement of the basilar membrane (BM), thus increase the loudness range (dynamic range) of the ear;
- increase frequency selectivity at the basilar membrane;
- protect inner hair cells by stiffening upon loud signals (until damaged).



**Békésy György**, 1899-1972  
Biophysicist  
Nobel-prize 1961 awardee


## •• Equal-loudness contours





- 1933: Fletcher & Munson (*JASA*)
- 1956: Robinson & Dadson (*BJAP*)
- 2003: Multisite study, averaged over many normal hearing (NH) listeners → **ISO 226:2003**
- At 1000 Hz: dB SPL  $\approx$  loudness.
- At other frequencies, e.g. at 50 Hz (mainly at higher levels) not.

Loudness perception is  
**frequency dependent**  
and **non-linear**.



## •• Scales and units in audio signal processing: dB HL



## •• Scales and units in audio signal processing: dB HL





## •• Scales and units in audio signal processing: dB HL



Bisgaard standard audiograms:

- S1: very mild (*degree of hearing loss*)



## •• Scales and units in audio signal processing: dB HL



Bisgaard standard audiograms:

- S1: very mild (*degree of hearing loss*)
- S2: mild

## •• Scales and units in audio signal processing: dB HL





Bisgaard standard audiograms:

- S1: very mild (*degree of hearing loss*)
- S2: mild
- S3: severe

Source: [Bisgaard et al.](#) (2010)

## •• Scales and units in audio signal processing: dB HL




Bisgaard standard audiograms:

- S1: very mild (*degree of hearing loss*)
- S2: mild
- S3: severe
- D: hearing loss bordering on deafness

Source: [Bisgaard et al. \(2010\)](#)




## •• Hearing and hearing impairments



## •• Hearing and hearing impairments

Questions:

- How many are affected?
  - [WHO 2021](#): worldwide ca. 466M people with disabling hearing loss.
- How bad is impaired hearing?
  - It is bad, especially in background noise. Try for yourself → see **slide** with hearing loss simulators.
- What are the causes of the impairments?
  - Genetic causes, complications at birth, infectious diseases, chronic ear infections, certain drugs, ageing, and exposure to excessive noise (“1.1B young people under 35 due to high SPLs in recreational settings”).
- What can be done to help people?
  - Hearing protection and education, access to early and on-demand screening;
  - access to hearing systems: hearing aids, cochlear implants, and other assistive devices ([WHO 2021](#): currently only 17% of those who could benefit from use of a hearing aid actually use one);
  - sign language and social support.



## •• Causes of hearing impairment



- Malformation
- Foreign body




- Eardrum scars / malformation / perforation
- Barotrauma (can be temporary)
- Otitis media or Eustachian tube dysfunction (fluid in the middle ear)
- Rigid joints between ossicles, otosclerosis



- Aging (presbycusis)
- Exposure to loud noise
- Trauma / Barotrauma / Head injury
- Genetic or hereditary hearing loss
- Viral infection of the hearing nerve
- Malformation of the inner ear or tumor
- Ototoxic drugs (antibiotics, some anti-inflammatory drugs, chemotherapy)



- Impacted ear wax
- Foreign body or fluid
- Inflammation or infection
- Malformation or tumor



## •• Causes of hearing impairment





- Malformation
- Foreign body



- Eardrum scars / malformation / perforation
- Barotrauma (can be temporary)

Exposure to high sound pressure level → **Hair cell damage**




- Foreign body or fluid
- Inflammation or infection
- Malformation or tumor

- Ototoxic drugs (antibiotics, some anti-inflammatory drugs, chemotherapy)

Source: [Ryan, 2000](#)

## •• Consequences of hearing impairment


- Communication is difficult to follow, uncertainties become frequent in everyday life.
- Environmental sounds missing (bird song, water dripping, water boiling over when cooking, ...).
- Withdrawal, due to
  - missing what others are talking about → frustration;
  - unwilling to ask repeatedly → shame, embarrassment;
  - everyday situations become challenging;
  - listening to speech becomes very effortful.
- Most of the time no sudden decline but a creeping process → unconscious adaptation.
- Often leads to "pro-active strategy": responses to what was meant to be said.
- May lead to loneliness, isolation, depression (false sense of others often being angry with them).
- **Unavoidable cognitive decline if untreated.** [FJ23]
- Danger: dementia. [FJ23]



[FJ23] F. Jiang et al., "Association between hearing aid use and all-cause and cause-specific dementia: an analysis of the UK Biobank cohort," *The Lancet*, vol. 8, no. 5, doi: 10.1016/S2468-2667(23)00048-8, 2023.


## •• First signs of a hearing impairment #1

- Speech at moderate level not always understood, but **loud speech becomes disturbing** and misinterpreted as offensive.



## •• First signs of a hearing impairment #1



- Speech at moderate level not always understood, but **loud speech becomes disturbing** and misinterpreted as offensive.



- In most cases, the amount of energy arriving at the IHCs is comparable to that of the NH (normal hearing) case.
- **Threshold of pain** is **unchanged** (and also the range of uncomfortably loud signals).
- The **dynamic range** (between threshold of hearing and threshold of pain) is **narrower**.
- In hearing assistive technology:
  - we need to **know the thresholds** (→ **pure tone audiometry**),
  - we need to **apply dynamic compression**.

## •• First signs of a hearing impairment #2

- Understanding **speech in noise** (other talkers, traveling, background music) becomes increasingly **difficult**.
- If pure tone audiogram looks OK → **hidden hearing loss**.
- Speech in noise test can clarify the issue.  
In German-speaking regions: OLSA (Oldenburger Satztest).



- Problems with understanding speech in noise often arise before a pure tone audiogram indicates hearing impairment.  
**Speech reception threshold, SRT, increases.**
- Threshold and equal loudness curves may remain largely unchanged.
- Frequency and amplitude resolution declines, which makes it the auditory system harder to distinguish signal from noise.
- Hearing assistive technologies may help with noise reduction algorithms.



## •• Hearing aids – Early history

Source: [hearingaidmuseum.com](http://hearingaidmuseum.com), 2021

Ear trumpets



17<sup>th</sup> century

Vactuphone  
(1<sup>st</sup> vacuum-tube HA)



1902

1921

Behind-the-ear  
hearing aids



1952

1956

Widex Senso  
full digital BTE HA



1984

1996



Hutchison's Acousticon  
(1<sup>st</sup> electrical HA)



Sonotone Model 1010  
(1<sup>st</sup> transistor HA)



Nicolet Phoenix  
(programmable digital HA)

## •• Hearing aids – Modern form factors



Behind the ear (BTE)



Receiver in canal (RIC)



In the canal (ITC)



Completely in canal (CIC)

Source: [betterhearing.org](https://betterhearing.org), 2021



- + least earwax impact
- definitely visible



- + least feedback



- + glasses-friendly
- ear occlusion



- + glasses-friendly
- ear occlusion
- usability

## •• Hearing aids – Requirements

- Attested hearing impairment is a **disability**.
- Health insurance (in Germany) pays for an adequate hearing aid system.

### **Technical requirements** (as of 2023):

- digital technology;
  - a minimum of 6 frequency channels with compression;
  - a minimum of 4 selectable, individualized programs (=settings);
  - directional microphones (to help separate speech and noise);
  - availability of feedback suppression method.
- 
- Plus, requirements of the **Medical Device Regulation (MDR)** of the European Union:
    - manufacturer must implement quality and risk management systems;
    - performance, user benefit, and user safety must be proven in clinical trials.

### **Disability:**

35+ dB hearing loss in  
better ear confirmed  
by a medical doctor

## •• Co-payments for hearing systems in Germany

Co-payment for hearing aid provision has existed in Germany since the health reform of 1989 (effective from 1993).

In 2023, the statutory health insurance funds [GK22] generally cover a contract price of

- €685 per hearing aid, plus
- a standard fee for custom-made earpieces of €33.50 and
- a flat-rate service fee for repair work of approximately €125.

Contract price generally covers standard hearing aid models (→ *Nulltarif*), but patients can pay extra to purchase premium products.

Cochlear implants are usually covered 100% by the health insurance.

[GK22] GKV-Hilfsmittelverzeichnis, "Festbeträge für Hörhilfen," 2022

## •• Hearing aids – Expected features

### General:

- Small, lightweight, intuitive, discrete\*.
- Remote control, smartphone app.
- Wireless streaming (TV, remote microphone).
- Quickly rechargeable, long-lasting battery.

### Audiological features:

- Low noise, high gain, high bandwidth (10 kHz+).
- 10+ frequency channels for DRC.
- Adaptive microphone directionality.
- Noise reduction: ambient, wind, impulse.
- Environment classification.
- Advanced feedback cancelation.
- Frequency shifting.
- Low latency ( $\ll 10$  ms), binaural algorithms.

Discrete

vs.

Fancy



## •• Hearing aids – Medical features

### General:

- Small, lightweight, intuitive, discrete.
- Remote control, smartphone app.
- Wireless streaming (TV, remote microphone).
- Quickly rechargeable, long-lasting battery.

### Audiological features:


- Low noise, high gain, high bandwidth (10 kHz+).
- 10+ frequency channels for DRC.
- Adaptive microphone directionality.
- Noise reduction: ambient, wind, impulse.
- Environment classification.
- Advanced feedback cancellation.
- Frequency shifting.
- Low latency ( $\ll 10$  ms), binaural algorithms.

### Fitness and medical features:


- Fitness tracking (step counter, pulse rate etc.).
- Blood pressure monitoring.
- On-demand activity log.
- Fall detection.



## •• Hearing aids – What is inside?




## •• Hearing aids – What is inside?







## •• Hearing aids – What is inside?




## •• Hearing aids – What is inside?




## •• Hearing aids – What is inside?



## •• Hearing aids – Basic signal flow



## •• Hearing aids – Architecture of a hearing aid ASIC



## •• Hearing aids – Basic specifications

**Rough estimates!**



Source: [jabra.com](https://jabra.com), 2021



Source: [amd.com](https://amd.com), 2021

	Hearing aid ASIC	Hearable ASIC	PC CPU
<b>Typical power consumption (algorithms on, wireless off)</b>	~ 1-2 mW	~ 15-50 mW	~ 35-125 W
Typical power consumption (algorithms on, BLE in use)	+ 3-8 mW	+ 15-50 mW	
<b>Clock speed</b>	~ 10-60 MHz	~ 20-200 MHz	
MIPS @ max. clock speed	~ 400 MIPS	~ 1 000 MIPS	~ 200 000 MIPS
<b>Typical dimensions</b>	~ 4 x 6.5 x 1 mm = 26 mm <sup>3</sup>	~ 5 x 8 x 1 mm = 40 mm <sup>3</sup>	~ 40 x 40 x 6 mm = 9.6 cm <sup>3</sup>
Number of cores	~ 2-8	~ 4-8	~ 2-32
Fabrication process	~ 45-22 nm	~ 45-22 nm	~ 14-7 nm


## •• Hearing aids – Construction of a hearing aid



Source: [kind.com](https://www.kind.com), 2021





Source: [taoglas.com](https://www.taoglas.com),  
2024




Source: [kind.com](https://www.kind.com), 2021

## •• Excursion: Telecoil




Source: [signet-ac.co.uk](http://signet-ac.co.uk), 2021






### T-coil or induction loop receiver:

- Operating at audio frequencies, no RF.
- Induction loop parts:
  - constant current amplifier,
  - closed loop cable.
- Attention (privacy):
  - EM-field bleeding.
  - Always broadcast.





## •• Next time





Thank you very much!  
Questions?

---

Dr.-Ing. Tamas Harczos  
[tamas.harczos@audifon.com](mailto:tamas.harczos@audifon.com)

audifon GmbH & Co. KG