



Hearing Systems – Part 3 of 3

<https://hearingsystems.github.io/>

TU Ilmenau – Audio Signal Processing & Audio Systems

January 25, 2023

Dr.-Ing. Tamas Harczos

•• Recap: Hearing aids – Features

General:


- Small, lightweight, intuitive, discrete.
- Remote control, smartphone app.
- Wireless streaming (TV, remote microphone).
- Quickly rechargeable, long lasting battery.

Audiological features:


- Low noise, high gain, high bandwidth (10 kHz+).
- 10+ frequency channels for DRC.
- Adaptive microphone directionality.
- Noise reduction: ambient, wind, impulse.
- Environment classification.
- Advanced feedback cancelation.
- Frequency shifting.
- Low latency (< 10 ms), binaural algorithms.

Fitness and medical features:


- Fitness tracking (step counter, pulse rate etc.).
- Blood pressure monitoring ([Sonion BiometRIC](#)).
- On-demand activity log.
- Fall detection.




•• Hearing aids – What is inside?



•• Hearing aids – What is inside?




•• Hearing aids – What is inside?




ASIC: application specific integrated circuit, HMI: human-machine interaction, MUX: multiplexer, NFMI: near-field magnetic induction, NVM: non-volatile memory, PM: power management, T-Coil: Telecoil (receiver for inductive loop), WL: (other) wireless technology

•• Hearing aids – What is inside?




ASIC: application specific integrated circuit, HMI: human-machine interaction, MUX: multiplexer, NFMI: near-field magnetic induction, NVM: non-volatile memory, PM: power management, T-Coil: Telecoil (receiver for inductive loop), WL: (other) wireless technology

•• Hearing aids – Basic signal flow



•• Hearing aids – Architecture of a hearing aid ASIC



ASIC: application specific IC, DAI: digital audio input, EEPROM: electrically erasable programmable read-only memory, GPIO: general purpose input/output, I²C: inter-integrated circuit, IOC: input/output controller, IP: intellectual property, LSAD: low-speed A/D converter, PCM: pulse code modulation, SPI: serial peripheral interface, UART: universal asynchronous receiver transmitter, WL: (other) wireless technology

•• Hearing aids – Basic specifications

Rough estimates!




Source: jabra.com, 2021



Source: amd.com, 2021


	Hearing aid ASIC	Hearable ASIC	PC CPU
Typical power consumption (algorithms on, wireless off)	~ 1-2 mW	~ 15-50 mW	~ 35-125 W
Typical power consumption (algorithms on, BLE in use)	+ 3-8 mW	+ 15-50 mW	
Clock speed	~ 10-60 MHz	~ 20-200 MHz	
MIPS @ max. clock speed	~ 400 MIPS	~ 1 000 MIPS	~ 200 000 MIPS
Typical dimensions	~ 4 x 7 x 1 mm = 28 mm ³	~ 5 x 8 x 1 mm = 40 mm ³	~ 40 x 40 x 6 mm = 9.6 cm ³
Number of cores	~ 2-8	~ 4-8	~ 2-32
Fabrication process	~ 45-22 nm	~ 45-22 nm	~ 14-7 nm

- Hearing aids – Construction of a hearing aid





Source: [kind.com](https://www.kind.com), 2021

•• Hearing aids – Construction of a hearing aid




Source: [kind.com](https://www.kind.com), 2021

•• Excursion: Telecoil




Source: signet-ac.co.uk, 2021





T-coil or induction loop receiver:


- Operating at audio frequencies, no RF.
- Induction loop parts:
 - constant current amplifier,
 - closed loop cable.
- Attention (privacy):
 - EM-field bleeding.
 - Always broadcast.



•• Recap: Equal-loudness contours




•• Hearing aids – Wide dynamic range compression



Linear amplification alone:

- Increases amplitude always with the same gain.
- Not conform with loudness perception of hearing impaired. (Will often cause uncomfortably loud percept.)
- Can lead to more hearing impairment.


•• Hearing aids – Wide dynamic range compression



Compression:

- Example:
40 dB increase in the input leads to 20 dB increase in the output, i.e.,
→ compression ratio: 40 dB : 20 dB = 2:1
- Compression ratio can vary, but usually does not exceed 4:1 in the hearing aid field.


•• Hearing aids – Wide dynamic range compression



Compression:

- Low level gain: gain at ca. 50 dB SPL input
- High level gain: gain at ca. 90 dB SPL input


•• Hearing aids – Wide dynamic range compression



Second linear range:


- Can be useful to allow the listener to realize that the input is very loud.
- Should be verified to not extend into the threshold of pain. (See “limiter” on next slide.)
- Remember the normal hearing ear handles loud noises very similarly! Basilar membrane nonlinearity:

Reminder



Adapted from: [Johnstone et al., 1986](#)


•• Hearing aids – Wide dynamic range compression



Limiter:

- Ensures not to cross the threshold of pain.
- Ensures technical limits not to be exceeded.
(MPO: maximum power output of a hearing device. Limited for example by the power source, amplifier, and receiver.)
- Is equivalent to a compressor with $\infty:1$ compression rate.


•• Hearing aids – Wide dynamic range compression



Expansion:


- Usually employed in the low SPL range.
→ Less gain added to soft noises.
- Reduces inherent noise of the system (e.g. that of microphones or preamps).
- Reduces irrelevant background noise.
→ Increases overall SNR.
- Typical expansion rate: 1:1.5 - 1:3.
(1:1.66 in the graph on the left.)

•• Hearing aids – Wide dynamic range compression




Gain plot:

- A shorter / compacter way to describe compression-expansion curves.




•• Hearing aids – Multi-channel WDRC



Adapted from: [Kiessling & Steffens, 1991](#)

Multi-channel compression:


- The normal dynamic range of speech varies across frequency.
- The amount of hearing loss also varies across frequency, individually.
- Solution: use multiple dynamic compressors, each fitted for a distinct frequency range.



- A computationally cheap, but high quality filter-bank with almost perfect reconstruction is needed.
- (Keywords for searching further information: FFT-based WOLA filter-banks and warped filters.)

•• Recap: Causes of hearing impairment

- Malformation
- Foreign body
- Eardrum scars / malformation / perforation
- Barotrauma (can be temporary)
- Otitis media or Eustachian tube dysfunction (fluid in the middle ear)
- Rigid joints between ossicles, otosclerosis




- Impacted ear wax
- Foreign body or fluid
- Inflammation or infection
- Malformation or tumor

- Aging (presbycusis)
- Exposure to loud noise
- Trauma / Barotrauma / Head injury
- Genetic or hereditary hearing loss
- Viral infection of the hearing nerve
- Malformation of the inner ear or tumor
- Ototoxic drugs (antibiotics, some anti-inflammatory drugs, chemotherapy)


Source: www.lobec.ca, 2021.

•• Cochlear implants




Source: [medel.com](https://www.medel.com), 2016

•• Cochlear implants




Source: [cochlear.com](https://www.cochlear.com), 2021

•• Cochlear implants



Source: [cochlear.com](https://www.cochlear.com), 2021 (unprocessed sound)


•• Cochlear implants – Signal processing strategy



Source: [G. Clark, Cochlear Implants](#), 2003

ACE signal processing strategy
(advanced combination encoder)

Biphasic, charge balanced pulses



Source: [Tognola et al.](#), 2005

•• Cochlear implants – Early history

Alessandro Volta
50V DC ear-to-ear



Lundberg (1950),
Djourno & Eyries (1957)

Direct electrical stimulation of the auditory nerve with distinguishable percepts in human.

House's first clinical
multi-channel CI



Source: [Clark, 2013](#)

Clarion implant
Advanced Bionics
(USA)



Source: [cochlear.com, 2021](#)

1790s

1868

1950-57

1961


1982

1985

1996


2001

Source: [abebooks.com, 2021](#)



Rudolf Brenner's
AC stimulation

Source: [Eisenberg, 2005](#)



William F. House's
single-channel CI

Nucleus implant
Cochlear Ltd. (AU)

Combi40+ implant
MED-EL GmbH (AT)

Commercialization →

•• Hearing aids (HA) vs. Cochlear implants (CI)

- System size

BTE processor is about the same size for both, but CI has more system parts (RF link, implanted stimulator, implanted electrode array).

- Power consumption

Signal processing requires about the same power for both, but CIs need extra current for electrical stimuli and RF link.

- Core signal processing

Both do spectral decomposition via an analysis filter-bank, but while HAs need to synthesize audio on the output, CIs calculate parameter for electrical stimuli.

- Additional algorithms / features

Very similar: directionality, dynamic range compression, noise reduction, binaural algorithms, but CIs do not need feedback canceller. Both can be streamed to wirelessly and be remote controlled.

- Current issues

Both: speech in noise enhancement, power consumption and size, connectivity, user-friendly handling.
CIs: reverb, **pitch discrimination** (12/16/22/24 channels, does not matter ← **current spread**).

•• Cochlear implants – Ways to improve pitch sensitivity

- Train cochlear implant users

Use the brain's plasticity. Needed anyway in CI therapy, but has its limits.

- Enhance pitch changes


Smart preprocessing can increase speech melody changes in real-time ([PREX project](#), [Fraunhofer IDMT](#)).



- Improve temporal pitch cues

Better stimulation strategies ([SAM project](#), [Dissertation](#), [Fraunhofer IDMT](#))

•• Pitch perception



•• Cochlear implants – Ways to improve pitch sensitivity

- Train cochlear implant users

Use the brain's plasticity. Needed anyway in CI therapy, but has its limits.

- Enhance pitch changes

Smart preprocessing can increase speech melody changes in real-time ([PREX project](#), [Fraunhofer IDMT](#)).



- Improve temporal pitch cues

Better stimulation strategies ([SAM project](#), [Dissertation](#), [Fraunhofer IDMT](#))




- Improve place pitch cues (reducing current spread)


- Better electrode design (“modiolus-hugging” or “adhering” electrodes),
- Current steering (parallel stimulation over multiple electrode contact points)
- Optogenetic stimulation (optical cochlear implants).

•• Optogenetic cochlear implants


Electrical stimulation



Optical stimulation



Research hardware




Challenges of optogenetic cochlea stimulation:

- Light-sensitive neurons: efficient, reliable and safe opsin* expression via *genetic manipulation* of the spiral ganglion neurons.
- Stimulation technology must be efficient, reliable and safe.

*opsin = light-gated ion channel

•• Optogenetics



Source: be.mit.edu, 2021

- Optogenetics – Demo with the research hardware



Music: Josh Armistead - Lazer Beams

Licence: Attribution 4.0 International (CC BY 4.0)

•• The Embedded Signal Processing (ESP) Team of audifon in Ilmenau

- Primary goals (among all) are:
 - monitoring hearing aid development trends;
 - preparation for future hearing aid HW platforms;
 - developing algorithms that may become part of the firmware of future hearing aids.
- Status quo:
 - Three full members, still missing an embedded C programmer.
 - Our first success: the *PrimeHA* project.
- PrimeHA:
 - Diploma work of Amelie Hintermaier entitled “*Entwurf einer kompakten, energieeffizienten, frei programmierbaren Hörgeräteplattform zur praxisnahen Erprobung audiologischer Algorithmen*”.
 - Got very good grade with distinction.



Tamas Harczos

Amelie Hintermaier

Iko Pieper

•• PrimeHA: a platform to develop and test hearing aid algorithms

Main features:

- ARM Cortex M4f @ 64 MHz, 256 kB RAM, 1 MB NVM.
- Current consumption ~3-10 mA with BLE 5.1 radio.
- Audio input: 2x PDM MEMS microphones.
- Extras: motion sensor and μ SD-card slot on board.
- Made of commercial off-the-shelf components.
- Fully C-programmable under free-for-all license.

Current use:

- Firmware framework development.
- Testing and optimizing audio signal processing algorithms like feedback canceller, multi-band compressor, noise reduction.
- Learning lessons on acoustically engineered mechanical design.





Thank you very much!
Questions?

Dr.-Ing. Tamas Harczos
tamas.harczos@audifon.com

audifon GmbH & Co. KG